

A high sensitivity, tumor-informed liquid biopsy platform, designed to detect minimal residual disease at part per million resolution

Sean Michael Boyle, Gabor Bartha, John Lyle, Jason Harris, Josette Northcott, Dan Norton, Rachel Marty Pyke, Fábio C. P. Navarro, Christian Haudenschild, Rose Santiago, Robin Li, Chris Nelson, Yelia Huo, Manju Chinappa, Qi Zhang, John West, Richard O. Chen
Personalis, Inc. | 1330 O'Brien Dr., Menlo Park, CA 94025

#5163

Contacts:
Sean.Boyle@personalis.com
Richard.Chen@personalis.com


Introduction

While circulating tumor-derived DNA (ctDNA) is an emerging biomarker for many cancers, the limited sensitivity of current detection methods reduces its utility for diagnosing MRD across a variety of clinical applications. Standard-of-care (SOC) radiological-based technologies, including CT, PET and MRI scans, also remain limited in their ability to detect residual disease during or after surgical or systemic therapy due to the minimum tumor volume required. NeXT Personal™, an advanced, personalized, and tumor-informed liquid biopsy assay, is designed to detect molecular residual disease (MRD) and cancer recurrence at the earliest timepoints — prior to, during or after treatment — in patients previously diagnosed with cancer.


Methods

NeXT Personal leverages tumor/normal whole genome sequencing (WGS) to design personalized, targeted MRD liquid biopsy panels for each patient. The MRD portion of the panel is composed of up to 1800 somatic tumor variants, enabling higher sensitivity MRD detection in plasma through tracking of high quality and lower noise variants. This allows the platform to achieve high sensitivity across cancer types and stages, including early-stage cancers and low mutational burden tumors, requiring only a single tube of blood (4mL plasma/15ng cfDNA), and 1mm³ of FFPE tumor tissue.


Results



References:
1. Minimum lesion detectability as a measure of PET system performance. Adler et al., EJNMMI Physics, 2017
2. Analytical validation of the Signatera™ MRD assay, a highly sensitive patient-specific multiple PCR NGS-based noninvasive cancer recurrence detection and therapy monitoring assay. Sethi et al., AACR, 2018, Abstract 4542
3. Genome-wide circulating tumor DNA monitoring for bladder cancer treatment management and organ preservation. Nordenskjöld et al., ASCO 2021, Abstract t5027
4. Phylogenetic tracking and minimal residual disease detection using ctDNA in early-stage NSCLC: A long TRACERx study. Akbani et al., AACR 2020, Abstract CT03
5. Analytical development of the RaBa® assay, a highly sensitive and specific assay for the monitoring of minimal residual disease. Maricic et al., AACR 2020, Poster 309
6. Phylogenetic ctDNA analysis depicts early-stage lung cancer evolution. Akbani et al., Nature, 2017
7. Detection of circulating tumor DNA in early- and late-stage human malignancies. Bettegowda et al., Science Translational Medicine, 2014



The MRD component of NeXT Personal leverages up to 1,800 WGS tumor derived somatic variants to provide a highly sensitive and aggregated measurement of tumor burden in plasma. Additionally, the ability to simultaneously track individual variants longitudinally can also be utilized to further our understanding of tumor biology and its dynamic response to therapy. NeXT Personal provides the ability to track and annotate individual variants over the whole evolution of a cancer patient's trajectory in a single panel design. Variants tracked by NeXT Personal in the blood are derived from those detected in the tissue, a Personalis-curated list of guideline-driven and resistance mutations, as well as those specified for panel inclusion by the user.



Central to the MRD performance of NeXT Personal is the number of tumor-derived mutations interrogated in patient plasma samples. Analytical and clinical MRD sensitivity can be increased by expanding coverage of loci containing tumor-specific mutations. Here, through an in-silico analysis, we highlight performance benefits that can be achieved by leveraging large numbers of high-quality somatic variants for MRD detection. Additionally, the platform performance is consistent using 15-50 ng of input material, which allows our platform to be applied with only a single tube of blood across a range of tumor types and stages.

NeXT Personal Yields Best-In-Class MRD Performance in Plasma, to a Lower Limit of 1-2 PPM




Figure 4: Cell line dilution series, HCC38, HCC1143, and HCC1937, demonstrates the linearity of tumor signal across dilution series down to the 1-3 PPM

We applied NeXT Personal to track MRD in a set of three well-characterized tumor-normal cell line systems. For each cell line, we serially diluted tumor cell line DNA with paired normal DNA, applying NeXT Personal to determine MRD signal (in PPM). Our platform reproducibly achieves an MRD limit of detection (LOD) of between 1 and 3 PPM. Furthermore, applying ddPCR, we corroborate our findings down to the detection limit of that technology.

NeXT Personal Effectively Tracks MRD Signal in Serially Diluted Patient Samples




Figure 5: NeXT Personal characterizes MRD LOD in a set of serially diluted patient samples

By serially diluting paired tumor-normal-plasma patient samples, we were able to reproducibly and systematically demonstrate NeXT Personal's ability to track variants down to the 1-3 PPM range on real patient samples. Here we demonstrate the relationship between dilution factor and detected PPM using both ovarian and an bladder cancer patient samples.

Applying NeXT Personal to Profile MRD For a Broad Set of Cancer Types and Stages




Figure 6: Patient T/N/P trios across 8 tumor types and 3 stages were sequenced demonstrating range of detection; 0.8- ~19,000 PPM

NeXT Personal Demonstrates 100% MRD Specificity When Tested on Healthy Donor Samples

Number of Healthy Donors	True Negative MRD Detections	False Positive MRD Detections	Specificity
40 Unique Individuals	40	0	100%

Figure 7: NeXT Personal detected no positive MRD detections when applied on healthy normal donors

We applied cancer patient panels to profile MRD in 40 healthy donor samples, NeXT Personal reported no false positive signals. Rather than using a predetermined minimum number of mutations to define MRD status, NeXT Personal uses advanced, proprietary statistical analyses to distinguish signal from noise for each individual patient panel and sample. Guided by a P-value threshold derived from a specificity target of >99.9%, NeXT Personal is optimized for both sensitivity and specificity.

Conclusion

NeXT Personal achieved highly sensitive and specific MRD detection, reproducibly demonstrating a LOD down to 1 PPM in different cancer types and cell line dilutions, representing approximately 10 to 100 times higher sensitivity than other liquid biopsy MRD approaches. The high sensitivity of NeXT Personal potentially enables MRD detection across a broad variety of cancers and stages, including typically challenging early stage, low mutational burden, and low-shedding cancers.

